



# Rocket<sup>®</sup> Colposcopy Range



**Rocket**medical

ISSUE 2



# Rocket® LLETZ Loop Range

The **Rocket®** LLETZ Loop Range has been designed to overcome the problems of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LLETZ electrodes.

- Unique 0.15mm wire versions to deliver higher cutting performance with significantly less thermal damage to the sample. (Small and Square styles only).
- 2.4mm (3/32") x 10cm shaft to fit all standard UK, US and European diathermy systems.
- Sterile - for single use. Guaranteed sterile and free of all contamination from previous use or reprocessing. (Unless the pack has been opened or damaged).

## A Rocket® Loops, round, sterile, single use

Small	14mm	25	R57800
Medium	18mm	25	R57801
Large	22mm	25	R57802
Square	10mm	25	R57803

## B DELTA™ Loops, half round, sterile, single use

Medium	18mm	25	R57851
Large	22mm	25	R57852

## C Rocket® NETZ™ Needles, sterile, single use

NETZ™ Needle	15mm	25	R57860
--------------	------	----	--------

## D Rocket® Coagulation Ball, round face, sterile, single use

Medium	5mm	25	R57804
Small	3mm	25	R57807

## Rocket® Diathermy Handswitch & Combo Packs

Diathermy Handswitch, sterile, 3m cable, US/UK fitting	1	R57865
ComboPack: Diathermy Handswitch, Rocket® Medium Loop, Rocket® 3mm coagulation ball	1	R57801-HS-03
ComboPack: Diathermy Handswitch, Rocket® Medium Loop, Rocket® 5mm coagulation ball	1	R57801-HS-05





**Rocket**medical



## rocketmedical.com

### UK Customer Services

Rocket Medical plc  
Sedling Road, Washington  
Tyne and Wear  
NE38 9BZ England  
Telephone: +44 (0) 191 419 6988  
Fax: +44 (0) 191 419 6989  
Email: [homesales@rocketmedical.com](mailto:homesales@rocketmedical.com)

### Head Office & Export Customer Services

Rocket Medical plc  
Imperial Way,  
Watford, Herts  
WD24 4XX England  
Telephone: +44 (0) 1923 651404  
Fax: +44 (0) 1923 240334  
Email: [exportsales@rocketmedical.com](mailto:exportsales@rocketmedical.com)

### US Office

Rocket Medical  
150 Recreation Park Drive, Unit 3  
Hingham, MA 02043, USA  
Customer Services: 1-800-707-ROCKET (7625)  
Telephone: 1-781-749-6223  
Fax: 1-781-749-6235  
Email: [usa@rocketmedical.com](mailto:usa@rocketmedical.com)