

Setting a New Standard for Recovery



Solace[®] System Features

- Advanced design that optimizes pain management
- Completely **latex-free** using patented synthetic rubber technology
- Optimized flow control with state-of-the-art *ZoneCatheter*[™]
- Compact and portable design including fanny pack
- No pump settings for patient or surgical staff
- *ZoneCatheter*[™] offered in multiple lengths

Solace[®] System Benefits

- The standard for performance and value in pain pumps
- Full surgical site coverage with *ZoneCatheter*[™]
- Peace of mind for patient and medical team
- Allows medical staff to focus on a patient's recovery – not their pain
- A better alternative than narcotics
- Improves recovery time and patient satisfaction



Targeted Coverage Ensures Consistent Delivery of Anesthesia

Not all catheters are the same.

The ZoneCatheter[™] (patent pending) is the result of extensive research and engineering. It incorporates a revolutionary new valve technology which allows for tighter spacing, versus standard holes in competitive catheters. With more valves per inch, ZoneCatheter[™] has a more consistent flow to ensure complete coverage of the entire surgical zone.

Solace[®] Specifications and Models

Specifications					Models				
Delivery Time (days)*	Nominal Fill Volume (ml)	Maximum Fill Volume (ml)	Nominal Flow Rate (ml/hr)	Number of Catheters	Standard Catheter	2.5" Zone Catheter [™]	5.0" Zone Catheter [™]	10.0" Zone Catheter [™]	Case Size (pkgs)
5 ½	65	75	0.5	single	S011	S021	–	–	5
2	100	125	2.0	single	–	S022	S032	–	5
5 ½	270	295	2.0	single	–	S023	–	–	5
3	270	295	4.0	dual	–	S024	S034	S044	5
4	400	500	4.0	dual	–	–	S036	S046	5
2	270	295	5.0	single	–	S025	S035	–	5
4	400	500	5.0	single	–	–	S037	S047	5
4	500	525	5.0	dual	–	–	S039	–	5
4	500	525	10.0	single	–	–	–	S048	5

* Delivery time is estimated based on nominal fill volume. Refer to Directions for Use for entire product information. Please contact Zone Medical for other Solace configurations.

Notes

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

THE SOLACE INFUSION SYSTEM IS STERILE WITH A NON-PYROGENIC FLUID PATH.

THE SOLACE INFUSION SYSTEM DOES NOT CONTAIN NATURAL RUBBER LATEX.

SINGLE PATIENT USE ONLY.

Distributed by Zone Medical, LLC



10064 Mesa Ridge Court, Suite 200
San Diego, California 92121
858-492-8570 Telephone
858-605-5919 Fax
888-259-0189 Orders

SolacePainRelief.com